



















Symbols Glossary

ISO 15223-1:2016 Medical Devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General Requirements		
Symbol	Symbol Title	Symbol Description
	Manufacturer	Indicates the medical device manufacturer.
	Date of Manufacture	Indicates the date when the medical device was manufactured.
	Lot Number (Batch Code)	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Serial Number	Indicates the manufacturer's serial number so that a specific medical device can be identified.
	Device Not Sterile	Indicates a medical device that has not been subjected to a sterilization process.
	Device Not Sterile	Indicates a medical device that has not been subjected to a sterilization process
	Sterilized with Irradiation	Indicates a medical device that has been sterilized using irradiation
	Do Not Re-Use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Use-by date	Indicates the date after which the medical device is not to be used
	Packaged Quantity	Indicates the amount of product included
	Material	Indicates material of the manufacturer
	Caution, Consult Accompanying Documents	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
 <small>ifu.spinalelements.com</small>	Electronic Instructions for Use	Consult Instructions for Use ifu.spinalelements.com
	Prescription Only	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.
	CE Marking	Complies with European Directives